



Cell Viability Incorporated

Reagent alternative for the Beckman Coulter Vi-Cell



Vi-Cell is a registered trademark of Beckman Coulter Inc.

[www.CellViability.com](http://www.CellViability.com)

# Proven Performance.....

A side-by-side comparison was done by a major pharmaceutical company with the Beckman Coulter reagents and the reagents from Cell Viability Inc..

Results proved to be comparable and acceptable.

Vi-CELL XR 2.03										
Beckman Coulter, Inc.										
Reagents Brand	Sample ID	Total cells	Viable cells	Viability	Total cells	Viable cells	Avg. diam.	Avg. circ.	Avg. cells	Avg. background
Beckman reagent	Molm13-1	1838	1821	99.1	2.03	2.01	13.85	0.86	36.8	204
	Molm13-2	1904	1881	98.8	2.10	2.08	13.52	0.87	38.1	204
	Molm13-3	2071	2050	99.0	2.29	2.26	14.17	0.88	41.4	204
Cell viability reagent	Molm13-1	1859	1843	99.1	2.05	2.04	13.39	0.89	37.2	205
	Molm13-2	1840	1821	99.0	2.03	2.01	13.22	0.89	36.8	205
	Molm13-3	2141	2115	98.8	2.36	2.34	13.74	0.89	42.8	206

**Conclusion: The reagents from the two brands have the similar results on cell counting with the same samples.**



# Part Numbers.....

Description	Cell Viability Inc Part Number	Beckman Direct Part Number	Comments
Cell Viability for Vi-Cell - Single Pack	2400102	383260	Beckman kit only brings 120 vials. Cell viability brings 200 vials. The reagent kit will run 200 tests. This is why Beckman users sometimes need to order extra vials
Cell Viability Kit for Vi-Cell XR - Quad Pack	2400101	383722	This kit from Beckman brings 8x120 vials to ensure enough vials for the 800 tests a Vi-Cell will get from a Quad Pak. Cell Viability kit brings 800 vials
Cell Viability Kit for Vi-Cell - Quad Pack	2400101	383198	This is an old part number that is still sold. Same reagents however the difference is that this kit brings 4x120 vials. Cell Viability kit brings 800 vials



# More Information.....

For more information, please contact Cell Viability Inc. at.....

[www.CellViability.com](http://www.CellViability.com)


[Sales@CellViability.com](mailto:Sales@CellViability.com)

1-800-618-7660

305-801-7140




**MORE THAN 50% OFF VI-CELL REAGENTS\***



The Cell Viability Reagent Kit from Cell Viability Incorporated is a direct replacement for Beckman Coulter Vi-Cell reagents at a significant savings. The performance of these kits have been proven to be of equivalent quality and performance by hundreds of satisfied customers.


**Part Numbers and Equivalents**  
**Cell Viability Single Reagent Pack** - p/n 2400169 includes 200 sample cups, equivalent to Beckman Coulter p/n 383260  
**Cell Viability Quad Reagent Pack** - a quantity of 4 of p/n 2400102 having a total of 800 sample cups, equivalent to Beckman Coulter p/n 383722 and 383198

For more information or ordering, please visit [www.CellViability.com](http://www.CellViability.com) or call us at 1-800-618-7660



**Cell Viability Incorporated**

8724 Sunset Drive, #403 Miami, Florida 33173  
Phone: 1-800-618-7660. Fax: 305-574-7800,  
e-mail: [Sales@CellViability.com](mailto:Sales@CellViability.com)



\* More than 50% off is based on current Beckman Coulter published list prices  
Copyright 2013 Cell Viability Inc. All rights reserved. Vi-Cell is a registered trademark of Beckman Coulter Incorporated.

Vi-Cell is a registered trademark of Beckman Coulter Inc.

[www.CellViability.com](http://www.CellViability.com)

